

User Guide

Device

Controller

v1.0

Table of contents

| | |
|---------------------------------------|----|
| 1. Device Controller | 4 |
| 1.1 Inscoper Device Controller | 4 |
| 1.2 Warning and cautions | 5 |
| 1.3 Specifications | 6 |
| 1.3.1 Operating specifications | 6 |
| 1.3.2 Input & outputs | 7 |
| 1.4 System requirements | 8 |
| 1.5 Installation | 9 |
| 1.6 Legal | 10 |
| 1.6.1 FCC/IC certification | 10 |
| 1.6.2 Copyright | 11 |
| 1.6.3 Disclaimer | 12 |

1. Device Controller

1.1 Inscoper Device Controller

The Inscoper system architecture decouples the two primary functions of traditional acquisition software: user interaction and device control. By separating these roles, the system ensures optimal performance, stability, and control across diverse microscopy setups.

- **User Interaction (Inscoper Imaging Software):** Configures the acquisition sequence, receives and displays the acquired images, and saves the experimental data.
- **Device Control (Inscoper Device Controller):** Communicates directly with the microscopy hardware to execute the acquisition sequence with precise hardware timing.

The Inscoper Device Controller acts as a dedicated hardware appliance and serves as the central communication hub between the workstation and the microscope's components.

Mermaid Diagram

By offloading real-time device management to the Device Controller, the Inscoper Imaging Software operates completely independently of hardware constraints. Consequently, regardless of the microscope's complexity, the user interface remains streamlined, intuitive, and focused on experimental requirements rather than underlying hardware limitations.

1.2 Warning and cautions

1. Always ensure that the Inscoper Device Controller is powered on before starting the computer. If uncertain, restart the computer.
2. The equipment can only be powered through Safety Extra Low Voltage that also complies with the limits of 6.3.1/6.3.2 of IEC 61010-1:2010.
3. Never use cables longer than 3 meters to connect devices (except for the Ethernet cable).
4. Using the Inscoper Device Controller in any manner not specified by Inscoper may compromise its safety protections.

1.3 Specifications


1.3.1 Operating specifications

| Parameter | Current Version | Next Versions | | |
|--|--|---|---|---|
| Device Controller Model | M | S | L | XL |
| Weight | 3.050 kg (6.724 lb.) | 1.150 kg (2.20 lb.) | 4.450 kg (8.82 lb.) | 4.650 kg (8.88 lb.) |
| Dimensions (LWH) | 230 * 230 * 230 mm (9.06 * 9.06 * 9.06 in.) | 120 * 200 * 65 mm (4.72 * 7.87 * 2.55 in.) | 420 * 434 * 84 mm (16.53 * 16.93 * 3.15 in.) | 420 * 434 * 84 mm (16.53 * 16.93 * 3.15 in.) |
| Power input | 24 VDC @ 2.5 A | 24 VDC @ 2.5 A | 24 VDC @ 3.75 A | 24 VDC @ 3.75 A |
| Operation temperature range | From 10 to 40 °C (from 50 to 104 °F) | From 10 to 40 °C (from 50 to 104 °F) | From 10 to 40 °C (from 50 to 104 °F) | From 10 to 40 °C (from 50 to 104 °F) |
| Altitude | Maximum 2000 m (6561.68 feet) | Maximum 2000 m (6561.68 feet) | Maximum 2000 m (6561.68 feet) | Maximum 2000 m (6561.68 feet) |
| Operating humidity range (non-condensing) | From 30 to 85 % | From 30 to 85 % | From 30 to 85 % | From 30 to 85 % |
| Storage temperature range | From 0 to 50 °C (from 32 to 122 °F) | From 0 to 50 °C (from 32 to 122 °F) | From 0 to 50 °C (from 32 to 122 °F) | From 0 to 50 °C (from 32 to 122 °F) |
| Storage humidity range (non-condensing) | From 30 to 85 % | From 30 to 85 % | From 30 to 85 % | From 30 to 85 % |

 **Warning**

Indoor use only!

1.3.2 Input & outputs



| | Current Version | Next Versions | | |
|-----------------------|--|---|---|--|
| TYPE | M | S | L | XL |
| ANALOG OUTPUTS | 2x (0-5V), 1x (0-12V) | 1 | 8 | 16 |
| | DAC resolution 12 bits, No Programmable Output range | DAC resolution 14 bits, Sample rate 180 MS/s, Output range ± 10 V, 0-5 V, +-5 V | DAC resolution 14 bits, Sample rate 180 MS/s, Output range ± 10 V, 0-5 V, +-5 V | DAC resolution 14 bits, Sample rate 180 MS/s, Output range ± 10 V, 0-5 V, +-5 V |
| ILDA | Via External Controller (MaxILDA) | - | - | 1 |
| ANALOG INPUTS | 2x (0-5V), 1x (0-12V) | - | - | 16 |
| | ADC resolution 12 bits, No Programmable Input range | - | - | ADC resolution 16 bits, Sample rate 1 MS/s, Input range ± 2.5 V, ± 5 V, ± 10 V, ± 12.5 V |
| I/O | 6 | 4 | 18 | 18 |
| SERIAL PORTS | 5 | 1 | 4 | 8 |
| USB host | 5 | 1 | 6 | 6 |
| Computer | Windows 7/10/11 | Windows 7/10/11 | Windows 7/10/11 | Windows 7/10/11 |
| | 1920 x 1080 px | 1920 x 1080 px | 1920 x 1080 px | 1920 x 1080 px |

1.4 System requirements

| | Minimum requirements |
|------------------|--|
| Operating System | Windows 10 / 11 64 bit |
| RAM | 32 GB |
| Hard disk drive | 1 TB NVMe SSD |
| Processor | Intel 14th Gen (minimum 8 cores) |
| Graphics Card | NVIDIA GeForce RTX 2000 Series |
| Screen | 1920 x 1080 resolution - Dual screens or a large 32-inch screen (2560 x 1440 (QHD)) highly recommended |

1.5 Installation

The Inscoper Imaging Software should only be installed by Inscoper staff or authorized representatives. The customer and/or user can participate in the installation process if an Inscoper representative provides explicit consent. In the event of modification to the microscopy system using Inscoper equipment or software without Inscoper's consent, the company declines all responsibility for any resulting consequences.

To request the installation of the Inscoper Imaging Software on a microscope, you must complete the following three steps:

1. Send a list of all devices connected to the microscope to contact@inscoper.com or via the dedicated form available at www.inscoper.com.
2. Provide the Inscoper team with three days of access to install the system.
3. Test the microscope with your team: acquire images of your research samples, and compare them with previous results.

1.6 Legal

1.6.1 FCC/IC certification

Any changes or modifications to this equipment not expressly approved by Inscoper may cause harmful interference and void the FCC authorization to operate this equipment.

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

This device must be professionally installed.

1.6.2 Copyright

The copyright in this document and the associated drawings are the property of INSCOPER and all rights are reserved. This document and the associated drawings are issued on the condition that they are not copied, reprinted, or reproduced, and that their contents are not disclosed, except in connection with the authorized use of the system.

The publication of information in this document does not imply freedom from any patent or proprietary right of Inscoper or any third party.

INSCOPER and the INSCOPER logo are trademarks of INSCOPER Company (INSCOPER SAS - 3771 boulevard des Alliés, 35510 Cesson-Sévigné - France). INSCOPER includes technology covered by the following patents:

- US Patent No. US10330911,
- EP Patent No. EP3123149,
- FR Patent No. FR3019324,

This product is updated periodically, and revisions will be incorporated into new editions of the user documentation.

1.6.3 Disclaimer

The information contained in this manual is provided on an “as is” basis, without any warranties, conditions or representations of any kind, whether express, implied, statutory or otherwise, including, but not limited to, any warranties of merchantability, non-infringement or fitness for a particular purpose.

In no event shall Inscoper be liable for any loss, or for any direct, indirect, special, incidental, consequential, or other damages, regardless of cause, whether arising in contract, tort, or otherwise in connection with the use of the information provided herein.